

MEDITECH

Solar Powered Refrigerator

Robust Construction, Technology, Quality, Reliability and Repeatability for extreme Tropic Conditions



- 12 V - 24 V DC operating voltage featuring shut-off at low voltage
- Dual Mode - Can work as Freezer and Refrigerator
- Eco-friendly cooling agent (R-134a)
- Maintenance Free Brushless DC Compressor
- Adjustable interior temperature
- Very low energy requirement
- Highly-efficient cooling
- Low Maintenance
- Double gasket to keep long cooling
- Easy to installation



Environmentally friendly Solar Refrigerator powered by Meditech Technology

ISO 9001, ISO1400, ISO13485, CE and WHO-GMP certified.

Meditech[®]



Bio Fresh provides the ultimate technology for all health-conscious food lovers. With Bio Fresh, fruit and vegetables, meat, fish and dairy products retain their healthy vitamins, delicate aroma and appetising appearance much longer than in a conventional fridge compartment.

Bio-fresh



Why Meditech Solar Refrigerator

Discover Technology, Quality, Reliability and Repeatability



Meditech Solar Refrigerator are rugged and designed for extreme tropic climatic conditions. Save on costs with Meditech high efficiency refrigerators and freezers . It has exceptionally low energy consumption, requiring smaller, less expensive power systems and low operating expense. High quality construction provides excellent reliability and long life.

Applications

Remote Homes	Unreliable grid	Micro-enterprises	Boats and Marine	cabins
Remote Stores	Ice Making	Beverage Vending	Eco-Resorts	Farms
Medical Clinics	Missionaries	Disaster Preparedness	Mobile Vendors	caravans

Meditech[®]

Know your Solar Refrigerator



Equipment meets electrical safety specification such as that of IEC (Class I)

Advance LED based PLC Based Controller



LED INDICATION

1. Line In
2. Power
3. Comp On
4. Heater On
5. Battery On
6. Battery Low
7. Temp High
8. Temp Low
9. Power Fail
10. Sensor Fail
11. System On












Scrolling LIVE Data logger on LED screen

1. Temperature
2. Incoming Voltage
3. Ambient Temperature
4. Time in hours of revolution chart
5. Current Date
6. Current Time
7. Battery Voltage

User Friendly Settings










1. Date
2. Time
3. High Temp Alarm
4. Low Temp Alarm
5. Hysteresis
6. Compressor Delay

LED Indicators

Sl.No	Photos	Function Indication
1.		This LED is on when there is a power supply to the machine from the main source
2		This LED is on when the Controller is in working condition.
3		This LED is on when the compressor is in working condition
4		This LED is on when the heater is in working condition
5		This LED is on when the battery is in working condition
6		This LED is on when the battery is discharged
7		This LED is on when the inside compartment temperature of the refrigerator is high above the set temp
8		This LED is on when the inside compartment temperature of the refrigerator is Low below the set temp
9		This LED is on when there is power failure.
10		This LED is on when there is sensor fail
11		This LED is on when all functions of the System is in working condition

Keypad Functions



Sl.No	Photos	Function Indication
1		The Menu button is pressed to enter into MENU option.
2		The Enter button is pressed to enter into “ENTER” option.
3		The Up Arrow button is pressed to go Upwards in the menu option.
4		The Down Arrow button is pressed to go Downwards in the menu option.
5		The data logger is pressed to get the data from the PLC to the computer
6		The MUTE button is pressed to MUTE any Alarms and sounds
7		The Door Open LED is activated along with the buzzer sound indicating the Door Open function
8		The ON button is pressed to switch on the refrigerator
9		The USB is provided to connect the Refrigerator to the computer system.

Description	MSR100	MSR150	MSR200	MSR300	MSR400	MSR500
Capacity in Litres	100	150	200	300	400	500
Dual Temperature Range Refrigerator and Freezer	Freezer upto to -20 Degree C (settable) Refrigerator from 2 to 8 degree C (settable)					
Power Supply	12-24 volts					
Controller and Display	PLC based LED Controller					
Temperature setting	(+/-) 0.1° C					
Temperature Uniformity	(+/-) 1° C					
Refrigerant	R134a	R134a	R134a	R134a	R134a	R134a
No. of Lids	1	1	1	1	2	2
Tank Material	Embossed Aluminium				PP	
Overall Dimension (mm)	568x560x845	760x560x845	980x560x845	1055x735x835	1440x620x835	1640x620x835
Net weight (kgs)	40	44	48	60	77	86
No of baskets	3	3	4	5	7	9
Wheels	Heavy Duty swivel caster wheels					