

Micro-controlled donor couch designed to provide a comfortable position for the donor



The Meditech Blood Donor Couch is a fully upholstered donor chair, designed to provide a comfortable position for the donor, and features variable positioning for either arm, as well as reclining and upright body positions.

Parameter	Value Description
Size	Length: 60-70 inches ; Width: 30-40 inches ; Head height 44-48 inches. Foot height -36-40 inches. Seat Height – 25-27 inches. Motor Power supply 220-250 volts AC
Self Weight	65 Kilograms
Weight Bearing	Ability to Accommodate upto 150 Kilograms
Positions	(1) Head Low (2) Feet High (3) Intermediate positions.
Operation Control	Remote Control
Hand rest	2 Numbers ( left and right)
Power Consumption	0.9 watts. Single phase. 230 Volts, 50/60 Hz
Castor wheels	4 Lockable castor wheels are provided.
Weight Capacity	200 Kilograms
Compliance	73/23/EEC(LVD) 89/336/EEC (EMC) 93/68/EEC
CE	This marking reflects with council directives on EMC 89/336/EEC, 92/31/EEC, 93/68/EEC.

1. All metal Frame Construction
2. Power slide control provides smooth shifting from **head-low, feet-high position** to any **intermediate position**
3. Accommodate assembly for operator ease and solid arm support
4. Ambar assembly for operator ease and solid arm support
5. 25 inch seat allows waist-high working level.
6. Adjustable arm rests swing out for easy seating as well as move up and down
7. Premium, UL listed polyurethane foam
8. Top quality, contract grade non flammable vinyl covering
9. Has springs, flexilators and layered foam for comfort
10. Posses a tray for positioning of blood donation equipment and supplies.
11. Compact built with a pull out step base
11. Non Bright blue Colour

Quality Standard: CE and WHO-GMP Certified