



The advantages of a blood donor station by providing a one piece lying surface with more flexibility in its functions. Besides your patient's requirements for a comfortable sitting position – it also provides a very relaxing lying position

Facility for blood collection from both sides. Provides a comfortable position for the donor. Variable positioning for either arm with comfortably wide arm-rests. Tilt adjustment can be done by using remote control. Comfortable chair type with soft padding for cushioning and rexin cover. Three Geared Actuator. Ensuring safety and comfort to the donor 4 Heavy Duty Lockable castors for easy mobility. Based on haemodynamic principles. (Donor's head can be lowered immediately and legs can be lifted above heart level so that blood can flow back to the brain and other vital organs, in case of vasovagal attack)

Catalogue No.	MTDC
Power Supply	230 VAC
Weight	80 KG
Lifting Capacity	200 kg
Control	Corded Remote
Power Consumption	100 W
Movement Actuation	DC Motor
Length of Head and Back Rest	80cm
Chairs total width from inside	60cm
Length of the seat	60cm
Length of the feet rest	50cm
Height of the chair base on the floor	55cm
Upholstery	Soft Upholstery of 4.0 or 10 cm thick
Actuator / Motors	3 nos