

ICE LINED REFRIGERATOR, MTBBR9



Special Features

- (1). Refrigeration System with Microprocessor controlled
- (2) The 670 Litres **ice lined refrigerator** is built with CFC-free polyethylene insulation.
- (3) Superior refrigeration from optimized evaporator and condenser design.
- (4) Dual interior layer design features *built-in ice packs* to keep chamber temperature within 2°C to 8°C range in an ambient temperature of no warmer than 43 °C for up to 24 hours during loss of power, presuming **ice-packs** are fully charged.
- (5). Operational ambient temperature range from 10°C~ 43 °C.
- (6) Temperature Control
- (7) The 670 Litres **ice lined refrigerator features** of digital temperature display.
- (8) Power on/off indicator.
- (9) Digital Temperature Display with accuracy of 0.1 °C
- (10) Inbuilt Battery back up of Temp display in case of power failure

Construction: Triple Walled

Inner Chamber: 304SS

Middle Chamber: Galvanised Steel Sheet

Outer Chamber: Powder coated GS Sheet.

Temperature Holding Time: 24 hrs

Insulation: PUF between inner and outer chamber minimum 75mm thick

Door: Solid Door, PUF Filled, Hinged at back with Lock, Latch and Key

Compressor: Hermetically Sealed Compressor

Refrigerant: Non-CFC, Non-HCFC

Sensors: Thermostat, Temperature Sensor with readout

Voltage: 230V 50Hz

Head Office: No.6 President Terrace. Takli Road, Dwarka. Nashik-422011 Phone:0253 - 6418632

Email: info@newmeditech.com. Website: www.newmeditech.com