

SPARES LIST



ICE LINED REFRIGERATOR

MTBBR

4. BILL OF MATERIAL

Serial Number	Part Description	Part Number
1	LED LAMP	MTBBRLL01
2	DOOR HINGES	MTBBRDH02
3	DOOR SWITCH	MTBBRDS03
4	SSR	MTBBRSSR04
5	CONTROLLER AND LED DISPLAY	MTBBRC05
6	ON/OFF SWITCH	MTBBROS06
7	PT100 SENSOR	MTBBRPTS07
8	COOLING COIL	MTBBRCC08
9	CONDENSOR FAN	MTBBRCF09
10	CASTOR WHEEL	MTBBRCW10
11	COMPRESSOR	MTBBRC11
12	DRYER UNIT DM15	MTBBRDU12
13	EXPANSION VALVE	MTBBREV13
14	AIR CIRCULATION FAN	MTBBRACF14
15	BASKET SS	MTBBRD15
16	DOOR HANDLE	MTBBRDH16
17	DOOR LOCK	MTBBRDL17