

## Euthanasia Chamber

Meditech Technologies India Pvt. Ltd is the largest manufacturer and exporter of Euthanasia Chamber. Euthanasia means that animals should be killed humanely with minimal pain and distress. Everyone has the responsibility to ensure that if an animal's life is to be taken, it is done with the highest degree of respect, and with an emphasis on making the death as painless and distress-free as possible. Taking these points into consideration Orchid has designed the Euthanasia Chamber

Meditech's euthanasia chamber includes the container, lid, Inlet valve with PTFE stopcock type controller, Outlet valve with PTFE stopcock type controller with adequate silicone tubing

Available different models suitable for small rodents and large rodents



### Features

- Simple quick and easy to use
- Bottom inlet decrease the induction time, "more humane"
- Complete Acrylic plastic cage – easy to see animals in chamber
- \* CO cylinder with regulator to be arranged by user

### Models

Model No.	Animal	Size (LXWXH)
EUTH 1	Mouse	7.5"x5.5"x4"
EUTH 2	Rat	11"x5.75"x7"
EUTH 3	Rabbit	16"x8"x10.5"
EUTH 4	Large animals	20"x12"x12"
EUTH 5	Large animals size 2	30"x12"x20"