

Fully Adjustable Acrylic Rabbit Restrainer

Meditech Technologies India Pvt. Ltd is the largest manufacturer and exporter of Fully Adjustable Acrylic Rabbit Restrainers. This restrainer is designed for use in two ways: In style 1, restrainer can be used with the help of curved back acrylic plate. This curved plate has a slider at bottom which can slide in the bottom sliders on restrainers according to animal size. In style 1, both height and length are adjusted using only one plate. In style 2, restrainer, height and length both can be adjusted using two plates. One plate as a cylindrical plate for height adjustment and second plate for length adjustment. Neck can be adjusted by sliding the front top plate and fixing it in the comfortable position neck slot. The floor of restrainer is specially made rough to avoid skidding.

Features

- Incorporates 2 designs in single instrument
- Ease of use and access to the rodent
- Made from thick perspex material.



Order Information

Model No.	Animal Weight Range	Overall Dimensions (mm)	Design Usage Style	Restraining Min Height (mm) from bottom	Restraining Max Height (mm) from bottom
LARIA	3 to 7.5 Kg	380 x 175 x 230	Style I (curved back plate)	220	470
			Style II (2 plate system- Curved top and straight back plate)	210	470