

## Guinea Pig Pans and Trolley

Meditech Technologies India Pvt. Ltd is the largest manufacturer and exporter of Guinea Pig Pans and Trolley.

### Guinea Pig Pans

Made of AISI non magnetic alkaline and chemical resistant 304 quality stainless steel material.

Pans are made of S.S 304 grade heavy duty 16 SWG sheet.

Pans are provided with 2 handles made from 4-5mm solid S.S. wire on each side for easy lifting

Pan has provision for holding 2 water bottles and supplied with 2 number of P.P. 500ml P.P. bottles with S.S. nozzle.

Each pan has holes of 20mm diameters on all sides

All the joints will be argon arc welded with smooth and glossy finish. Overall inner size: 920mm x 600mm x 230mm (LXWXH)



### Trolley/ Rack For Guinea Pig Pan

Made from 25mm x25mm x 16swg S.S. 304 square pipe

Suitable to hold 3 guinea pig pans with one pan in each shelf.

Movable on 4 No. of Polyurethane swivel castor wheels 2 with brakes.

The entire rack will be provided with stoppers on the backside of the angles to prevent cage falling from the back.

All the joints will be welded with argon arc welding with glossy finish

### Guinea Pig Pans - Ordering Information

Product	Model	Accessories Supplied
Guinea Pig Pan S.S.	GPP-1	S.S. 304 Pan with 500ml bottle and bottle holder

### Trolley / Rack To Hold Guinea Pig Cages - Ordering Information

Product	Model	Number of Pans
Rack for Guinea Pig Pan	GPP-T-1	3