

Trolley with Inbuilt Polypropylene Guinea Pig Cage

Meditech Technologies India Pvt. Ltd is the largest manufacturer and exporter of Trolley with Inbuilt Polypropylene Guinea Pig Cages. System has six guinea pig cages and each guinea pig cage comprises of P.P. tray forming bottom tray, S.S. door, top & side mesh inbuilt fitted on the trolley.

Polypropylene Tray: Tray is made of autoclavable food grade P.P. material and has provision for fixing food hopper. Size of P.P. tray is 650 x 500 x 200mm.

Rack : The rack structure will be fabricated from AISI 304 quality stainless steel square pipe of Size 20mm x 20mm. Sides of racks comprising cage top & side parts are made of 2.5mm diameter S.S. 304 rod.

Door : The door is fitted on the main frame of trolley and its size is 400mm x 195mm.

Food hopper : Fabricated from S.S. 304 sheet of size 130 x 100 x 145mm and fits on the polypropylene body.

Water bottle : Made of autoclavable food grade P.P. material of capacity 500ml with S.S. nozzle

Wheels : The entire battery with cages will move on 4 nos. of Polyurethane swivel 4” castor wheels two with brakes and two without brakes.

All the joints will be argon arc welded with smooth and glossy finish.
Overall Size: 1402 x 525 x 1610 (L x B x H) without wheels



Order Information

Product	Model	Accessories Supplied
Trolley with inbuilt P.P. guinea pig cage	GPINT-PP-01	Trolley with sides, top and door fitted,G.P. Cage tray- 6, each cage supplied with food hopper and water bottles