

## Water Feeding Bottles

Meditech Technologies India Pvt. Ltd is the largest manufacturer and exporter of Water Feeding Bottles.

For feeding water to experimental animals

Meditech Technologies produces a wide range of water bottles. The bottles are made of non transparent Poly-propylene or transparent Polycarbonate (PC).



### Features

- Leakage proof- No water leakage from bottle
- Assured water delivery- Each bottle cap is assured for proper water delivery
- Consistent thickness
- Wide mouth to aid filling and emptying, whilst also ensuring effective washing
- Volume graduations to aid volume measurements

## Polypropylene Water Bottle (Vacuum Style, Without Ball)

Press Fit Style: Cap made of autoclavable silicon bush with S.S. 304 protection flap

Threading Style: Cap is made of P.P. threading type with S.S. nozzle moulded in cap

### Models:

Model No	Capacity	Animals	Bottle	Cap	Sipper tube type
SI15P	150ml	Mice	P.P.	Silicon bush with S.S. Cap Press fit Style	S.S. 304 Nozzle
PT15P	150ml	Mice	P.P.	P.P. threading style	S.S. 304 Nozzle
SI25P	250ml	Mice, rat	P.P.	Silicon bush with S.S. Cap Press fit Style	S.S. 304 Nozzle
PT25P	250ml	Mice, rat	P.P.	P.P. threading style	S.S. 304 Nozzle
SI30P	650ml	Rat, Rabbit, Guinea Pig	P.P.	Silicon bush with S.S. Cap Press fit Style	S.S. 304 Nozzle
PT30P	650ml	Rat, Rabbit, Guinea Pig	P.P.	P.P. threading style	S.S. 304 Nozzle

## Polycarbonate Water Bottle (Vacuum style, without ball)

Press Fit Style: Cap made of autoclavable silicon bush with S.S. 304 protection flap

Threading Style: Cap is made of P.P. threading type with S.S. nozzle moulded in cap

### Models:

Model No	Capacity	Animals	Bottle	Cap	Sipper tube type
SI15C	150ml	Mice	P.C.	Silicon bush with S.S. Cap Press fit Style	S.S. 304 Nozzle
PT15C	150ml	Mice	P.C.	P.P. threading style	S.S. 304 Nozzle
SI25C	250ml	Mice, rat	P.C.	Silicon bush with S.S. Cap Press fit Style	S.S. 304 Nozzle
PT25C	250ml	Mice, rat	P.C.	P.P. threading style	S.S. 304 Nozzle
SI30C	550ml	Rat, Rabbit, Guinea Pig	P.C.	Silicon bush with S.S. Cap Press fit Style	S.S. 304 Nozzle