

## Algometer

Meditech Technologies India Pvt. Ltd is the largest manufacturer and exporter of Algometers.

Meditech Technologies's Algometers easily measures pressure thresholds and pain tolerances reported by patients. Algometry is ideal for helping evaluate fibromyalgia and trigger point tenderness.

Algometers are designed to quantify and document levels of tenderness via pressure threshold measurement and pain sensitivity via pain tolerance measurement. Pressure algometry is a reliable measure of pain in muscle, joints, tendons, and ligaments.



## Features

- Digital and Dial type (analogue) models available
- High precision and resolution
- Algometer PC Connectivity and software for analysis in case of digital models
- Digital Timer: 999 Min
- Supplied with tip of 1cm sq. diameter
- Set zero key
- In built rechargeable battery for digital model

## Models

Model	Dial type	Digital
Specification	ALGO-AN-01	ALGO-D-01
Capacity	20Kg	20Kg
Units	Kg, N	Kg, N, lbs
Load Division Value	200gm	0.01Kg
Accuracy	±0.5%	±0.5%
PC Connectivity and Software	No	Yes



## Order Information

Model	Certifications	Power requirements	Accessories
ALGO-AN-01	CE Compliant	N.A.	Dial type sensor with 1 tip of 1cm sq., carrying case
ALGO-D-01	CE Compliant	0/230V AC 50Hz 110/110V AC 50-60Hz*	Digital type sensor with 1 tip of 1cm sq., carrying case RS232 Cod, CD for software