

All Purpose Equipment

Meditech's All-purpose system is based on a powerful main drive unit. The various attachments can be easily and quickly connected by means & a universal gear box or directly. The AP-01 Plus is equipped with a single phase PMDC motor. The sturdy ball bearing worm gear drive is completely maintenance free. The unit features electronic speed control variable between 10 RPM to 300 RPM, A digital rpm controller regulates the speed. The All - Purpose equipment is ideal for small scale production in the pharmaceutical, chemical, cosmetic & food industries. It fits perfectly in to research & developments of tablets, ointments, creams & more. The modular design is compact, economical & offers a powerful motor drive with a wide range of easily interchangeable attachments for various purposes.



The different possibilities are as follows:

- Polishing
- Grinding
- Homogenizing
- Stirring
- Mixing
- Granulating
- Coating

Features

- Digital RPM Indicator: 10 - 300 RPM
- Digital Temp. Controller for Coating pan (5 °C above ambient to 150 °C)
- Emergency Switch
- Digital Timer: 999 Min
- Hot air gun arrangement for coating pan
- Auto and Manual Mode

Accessories

- Universal Gear Box
- Ball Mill
- Coating Pan
- Polishing Pan
- Double Cone Blender
- Cube Mixer
- V- Blender
- Homogenizer
- Pelletizer
- Planetary Mixer
- Dry Granulator
- Triple Roller Mill

