

All Purpose Equipment

Meditech Technologies India Pvt. Ltd is the largest manufacturer and exporter of Digital Enhanced Thin Layer Chromatograph de-TLC.


de-TLC is an inexpensive system, specially designed for PC-based documentation and analysis of thin layer chromatography (TLC) plates.

Scientist on a computer, using our exclusive de-TLC software, can quickly produce multi-spectral scans and calibration curves—output previously available only from more expensive equipment or complex procedures.

de-TLC combines an exceptionally versatile software module with a robust and easy to use imaging system to give you an outstanding instrument for all TLC plate applications.

Built to the highest standards but at an economic cost, de-TLC is the ideal choice for your laboratory.

Get de-TLC software and spend time on real research work instead of tedious measurements



Features

- Password protected system
- Accurately quantify the compounds
- Multiple image comparisons for statistical reliance
- Interactive graphical analysis
- Image refining
- Image documentation
- Keyboard shortcuts
- Interact with results and find relevant information within seconds instead of days.
- Inbuilt Calibration mode for periodic calibration
- Print and export the data
- Email the data for publications directly
- Create customizable reports
- Search stored reports