

## Isolated Tissue Bath

Meditech Technologies India Pvt. Ltd is the largest manufacturer and exporter of Isolated Tissue Bath.

Isolated organ bath are used to study effect of drug/chemical substances on isolated tissues in vitro.

Orchid's isolated tissue/organ bath systems are designed for in vitro investigations on isolated smooth, cardiac or skeletal muscle preparations.

Orchid's isolated organ bath systems are available as 1, 2 or 4 channel system and used in pharmacological research worldwide.



### Features

- No leakage- **Thick Perspex joint free tank**
- No glass parts- Organ bath tube is made of P. P. And reservoir coil is made of Perspex making system breakage free and maintenance free
- No leakage from organ tube fitting- Use of polypropylene bush for fixing organ tube in place of rubber cork
- Uniform Temperature- Inbuilt water circulation pump for uniform temperature control
- Available as 1, 2 or 4 channel system
- Available with analogue and digital temperature controller

## System Specifications & Models

Specifications	SE-1AW	SE-1D	SE-2A	SE-2D
Water Bath	One Bath	One Bath	One Bath	One Bath
Water stirrer	Yes	Yes	Yes	Yes
Temperature controller	Analogue	Digital	Analogue	Digital
Organ bath tube	01	01	02	02
Temperature range		5°C above ambient TO 60°C		
Temperature accuracy	±5°C	±5°C	±0.2°C	±0.2°C
Copler Set	01 set	01 set	02 set	02 set
Aeration tube	01	01	02	02
Power requirements	220/230V AC 50Hz 110/110v AC 50-60Hz*			

## Order Information

SE-1AW	SE-1D	SE-2A	SE-2D
Single unit with analogue temperature controller	Single unit with digital temperature controller	Double unit with analogue temperature controller	Double unit with digital temperature controller
Organ bath tube:1 Reservoir coil:1 Copler set:1 Aeration tube:1 Adequate Silicon tube for interconnections	Organ bath tube:1 Reservoir coil:1 Copler set:1 Aeration tube:1 Adequate Silicon tube for interconnections	Organ bath tube:2 Reservoir coil:2 Copler set:2 Aeration tube:2 Adequate Silicon tube for interconnections	Organ bath tube:2 Reservoir coil:2 Copler set:2 Aeration tube:2 Adequate Silicon tube for interconnections
Optional Accessories	Aerator,Tubing for Aerator,Marriot bottle 1 ltr,Stand for marriot bottle		