

## Spirometer

Meditech Technologies India Pvt. Ltd is the largest manufacturer and exporter of Spirometer.

A spirometer measures lung capacity and respiratory volume, it uses these measurements to create a spiograph. A spiograph represents the amount of air that moves into and out of the lungs with each we take. The spirometer records the volume and capacity of the lungs.

### Features

- Tank made of Stainless steel
- 6 liter capacity for vital capacity
- Use as a gas meter with chain compensated counter balance to the float balanced through cut its movement.
- Pulley calibrated to denote volume
- Side cock for gases, cock for draining off water
- Complete with corrugated tube mouth piece & recording lever.

