

BLOOD DONOR COUCH, MTDC04 With Foot Rest - Hydraulic System



The advantages of a blood donor station by providing a one piece lying surface with more flexibility in its functions. The Chair is equipped with Non Moving Parts and the whole position is changed using hydraulic system .Besides your patient's requirements for a comfortable sitting position – it also provides a very relaxing lying position. Designed for Comfort of the donor in Semi Sitting Position with Foot Support. Handles are provided to change the position of the chair for leg movement and Back Rest Movement.

Facility for blood collection from both sides. Provides a comfortable position for the donor. Variable positioning for either arm with comfortably wide arm-rests. Tilt adjustment can be done . Comfortable chair type with soft padding for cushioning and rexin cover. Ensuring safety and comfort to the donor 4 Heavy Duty Lockable castors for easy mobility. Based on haemodynamic principles. (Donor's head can be lowered immediately and legs can be lifted above heart level so that blood can flow back to the brain and other vital organs, in case of vasovagal attack)

Catalogue No.	MTDC04
Power Supply	NA
Weight of the chair	60 KG
Lifting Capacity of the chair	200 kg
Control	Hydraulic System
Movement Actuation	Hydraulic Lever System
Foot Support	Foot rest is provided additionally
Length	200 cm full chair in sitting position
Width	80cm
Upholstery	Soft Upholstery of 4.0 or 10 cm thick